

Turning Unsearchable Mixtures into Actionable Intelligence for Reliable DNA Evidence

**American Academy of Forensic Sciences
Criminalistics Section**

February, 2026
New Orleans, LA

**Edward Schikel¹, MBA; Terrance Lewis¹, MA;
Matthew Legler², BS; Mark Perlin², PhD, MD, PhD**

¹Montgomery County Detective Bureau, Pennsylvania
²Cybergenetics Corp, Pennsylvania



Cybergenetics © 2003-2026



Disclosure

- The views and opinions expressed in this presentation are those of the presenters & authors and do not represent any official views or opinions of the American Academy of Forensic Sciences.
- Some of the authors are employed by Cybergenetics.
- We show the TrueAllele® evidence interpretation technology developed by Cybergenetics.

HIGH VOLUME CRIMES
BURGLARIES, CAR BREAK-INS & GUN POSSESSIONS

- Burglary
- Thefts
- Car Break-Ins
- Weapons Possession

BURGLARIES CAR BREAK-INS

DNA Investigative Process

- BodeHITS Process
 - Enables DNA use in high-volume, low-level crime
 - Batch testing provides speed and cost efficiency
 - Regional database sharing links profiles across jurisdictions
-
- Trade-off
 - High-volume low-level crimes result in samples that are deemed non-searchable mixtures (NSM)



RapidHit Investigative Process

- Enables the 90-minute processing of single-source samples, producing DNA profiles for uploading and searching of the BodeHITS Database for regional sharing matches
- It allows us to answer investigative questions in real-time



Investigative Hurdle

Searchable vs. Non-searchable

- ✓ Searchable Profiles = Match or No match
- ✓ Non-searchable Mixtures = Stalled investigation

A New Question was asked

- Is there more we can do with the NSM Data on a monthly basis?



MontCo and Cybergenetics Collaboration



TrueAllele® Technology

Separates data into contributor genotypes
 uses all data, no thresholds, 10 contributors
 Represents uncertainty using probability
accurate: preserves identification information
unbiased: doesn't know suspect profile

Reports likelihood ratio match statistic
 change in genotype probability, error rates
 Automated computer workflow frees people
 no calibration, learns from evidence data

Number of Contributors



TrueAllele Database

Genotype matching database
 uploads all genotypes
 automated match server

Connects separated genotypes
 evidence to evidence
 evidence to reference

LR connection strength
 records DNA matches

TrueAllele + LIMS = TAID

TrueAllele Identification Database (TAID)
 uploads all contributor genotypes
 preserves DNA match information
 records all strong LR connections
 annotates with LIMS information


- Reduces human effort ten-fold
- Passes cost saving on to investigator
- Enables affordable volume crime DNA

Evidence to Evidence

	B						
Collected	14-Nov-22						
Door Handle	1						
Gear Shifter	15						
SteeringWheel	23						
log(LR) range	4-10	4-11	4-9	5-11	4-11	4-9	5-7
# Matches	9	39	7	4	9	6	4
	A	C	D	E	F	G	
Collected	31-May-21	6-Jan-23	13-Feb-23	16-Feb-23	22-Mar-23	2-Oct-23	
Door Handle		3					
Gear Shifter	3	2	4	3	4	2	
Ski Mask	3						
SteeringWheel	3	2		3	2	2	
Tire Iron				3			

Evidence to Reference

	R						
	Collected	10-Jan-23					
	Suspect	18 Y/O Male					
	log(LR) range	4 - 30					
	# Matches	18					
log(LR) range	23 - 27	4 - 9	17 - 30	18 - 30	16 - 30	18 - 20	16 - 21
# Matches	3	3	3	2	3	2	2
	A	B	C	D	E	F	G
Collected	31-May-21	14-Nov-22	6-Jan-23	13-Feb-23	16-Feb-23	22-Mar-23	2-Oct-23
Door Handle			1				
Gear Shifter	1	1	1	2	1	1	1
Ski Mask	1						
Steering Wheel	1	2	1		1	1	1
Tire Iron					1		



From Mixtures to Clarity

Montgomery Bode data to Cybergenetics

100 DNA items / month

60% evidence, 40% reference

⇒ 60% *unsearchable mixtures* unused

TrueAllele ID information improvement


1.32 evidence matches / item

98% of evidence is informative

⇒ 72% *unsearchable mixtures* matched


Case Example

- Juvenile male "HB" criminal trajectory
HB was the suspect in numerous vehicle thefts/
break-ins, but wasn't linked because of NSMs.
- HB was the suspect in numerous vehicle
thefts/break-ins, but wasn't linked because
of NSMs.
- A father and his two children were carjacked
by two suspects. One was masked, and the
other was identified and arrested.
- The two years of DNA were evaluated by
TAID, resulting in matches in seven cases
identifying HB, one of whom was the
masked perpetrator in the carjacking.



Montgomery County TAID

- Local PDs and Prosecutors have seen the revival
of stalled cases
- Associations are given on a monthly basis,
between cases in the same jurisdiction and
across jurisdictions
 - Evidence-Evidence
 - Evidence-Reference
- We realized that NSMs don't mean non-useable



Montgomery County TAID

Where do we go from here?

- Continue to prioritize DNA recovery
and analysis on volume crimes
- Support intervention before
escalation occurs
- DNA shouldn't just close cases—
It should change trajectories
of early offenders



More TAID information

<http://www.cybgen.com/information>



- Courses
- Newsletters
- Newsroom
- Presentations
- Publications
- Reliability
- Webinars

<http://www.youtube.com/user/TrueAllele>
TrueAllele YouTube channel